


About IPU (Image Processing Utilities) Plug-in

2.12

Introduction

Image Processing Utilities (IPU) is a plug-in for Gatan's DigitalMicrograph™.

This plug-in will provide compatible commands of Process menu regarding FFT as shown below:

With SHIFT key	With CTRL key
 The menu shows standard 2D processing commands: FFT (Ctrl+Alt+F), Inverse FFT (Ctrl+Alt+I), Reduced FFT, Live, Statistics, Live Statistics, Math, Complex To Real, Auto Correlation, Cross Correlation..., Convolution..., Wiener DeConvolution..., Extract Slice(s), Create Stacked Image from, Remove Spike Noise from SI..., Display Tags (Ctrl+Alt+D), About IPU..., and Help.	 The menu shows 3D processing commands: 3D FFT (Ctrl+Alt+F), 3D Inverse FFT (Ctrl+Alt+I), 3D Reduced FFT, Live, 3D Statistics, 3D Live Statistics, 3D Math, 3D Complex To Real, 3D Auto Correlation, 3D Cross Correlation..., 3D Convolution..., Wiener DeConvolution..., Extract Slice(s), Create Stacked Image from, Remove Spike Noise from SI..., Display Tags (Ctrl+Alt+D), About IPU..., and Help.

IPU Menu

Most of the commands switch to 3D processing by pressing the SHIFT key when the menu is selected, and to 2D processing for each slice by pressing the CTRL key.

These commands will work on an image with any image size, namely without a restriction of power of two, based on Intel's MKL (Math Kernel Library). We also provide with this plug-in the FFT-related functions that you can use from your script. These functions are also used internally by our several high-performance plug-ins.

This is a free plug-in, and we are happy to get your comments or suggestions on this plug-in.

Contents

-
- [Contact Us](#)
 - [Copyright Statements](#)
 - [Installation](#)
 - [Quick Reference](#)
 - [Functions](#)
-

Contact Us

General enquiries on the HREM Mouse should be sent to:

HREM Research Inc.
14-48 Matsukasedai
Higashimatsuyama
Saitama 355-0055
email: support@hremresearch.com
Website: <http://www.hremresearch.com/>

Generated on Thu Mar 29 16:48:24 2018 for Image Processing Utilities by  1.6.1

Copyright Statements

IPU

© Copyright 2008-2018 HREM Research Inc.

All rights reserved. This manual is protected by international copyright laws and treaties. Unauthorized reproduction and distribution of this manual, or any portion of it, will be prosecuted to the maximum extent possible and may result in severe civil and criminal penalties.

DigitalMicrograph is a trade mark of Gatan Inc.

Acknowledgements

We are grateful to Gatan software term to support our plug-in development.

Generated on Thu Mar 29 16:48:24 2018 for Image Processing Utilities by



Installation

1. Exit (Quit) DigitalMicrograph if it is launched.
 2. Copy IPU.dll and IPU.gtk to DigitalMicrograph/PlugIns.
-

Generated on Thu Mar 29 16:48:24 2018 for Image Processing Utilities by



Quick Reference

IPU Menu

		With SHIFT key	With CTRL key
FFT	Ctrl+Alt+F	3D FFT	Ctrl+Alt+F
Inverse FFT	Ctrl+Alt+I	3D Inverse FFT	Ctrl+Alt+I
Reduced FFT		3D Reduced FFT	
Live	▶	Live	▶
Statistics	▶	3D Statistics	▶
Live Statistics	▶	Live Statistics	▶
Math	▶	3D Math	▶
Complex To Real	▶	3D Complex To Real	▶
Auto Correlation		3D Auto Correlation	
Cross Correlation...		3D Cross Correlation...	
Convolution...		3D Convolution...	
Wiener DeConvolution...		Wiener DeConvolution...	
Extract Slice(s)	▶	Extract Slice(s)	▶
Create Stacked Image from	▶	Create Stacked Image from	▶
Remove Spike Noise from SI...		Remove Spike Noise from SI...	
Display Tags	Ctrl+Alt+D	Display Tags	Ctrl+Alt+D
About IPU...		About IPU...	
Help	▶	Help	▶

IPU Menu

Most of the commands switch to 3D processing by pressing the SHIFT key when the menu is selected, and to 2D processing for each slice by pressing the CTRL key.

FFT / 3D FFT

Same as DigitalMicrograph

Inverse FFT / 3D Inverse FFT

Same as DigitalMicrograph

Reduced FFT / 3D Reduced FFT

Same as DigitalMicrograph

Live

Same as DigitalMicrograph



Live Menu

FFT


Reduced FFT

Auto Correlation / 3D Auto Correlation

Same as DigitalMicrograph

Cross Correlation / 3D Cross Correlation

Extended cross-correlation function. By default this behaves same as DigitalMicrograph. However, if you choose this command with ALT key down, you will get a dialog shown below where you can select a cross-correlation mode and a mask size.



Cross Correlation Dialog

Mode

Cross correlation mode
1..Conventional, 2..Square Root Moduli, 3..Phase

Mask size


Mask radius given in pixels in Fourier space, or 0 to 0.5 in terms of Nyquist frequency.
Circular mask is used to remove high frequency noise.

Convolution / 3D Convolution

Calculates a convoluted image of selected two images.

Statistics / 3D Statistics

Same as DigitalMicrograph



Statistics Menu

Live Statistics

Use this command after placing a Rectangular ROI on an image. The selected statistics within the ROI will be dynamically updated when you move the ROI.

Sum
Mean and Std
Min and Max
All

Live Statistics Menu

Math / 3D Math

Same as DigitalMicrograph

a + b
a - b
a x b
a / b

Math Menu

Complex To Real / 3D Complex To Real

Creates a real image of the specified type from a complex image.

Phase
Amplitude
Real
Imaginary
Intensity

Complex To Real Menu

Extract Slice(s)

Extracts one currently displayed slice, or all slices, or selected slices from a front image.

Current Slice
All Slices
Selected Slices...

Extract Slice(s) Menu

Create Stacked Image


Creates a stacked image from the top most opened n-images in the image list of the window menu, or from the saved image in the selected folder, or from an open stacked image.

Opened Image...
Saved Images...
Stacked Image...

Create Stacked Image Menu

Remove Spike Noise from SI

Remove a spike noise(s) from a SI data. A spike noise is defined as the data point that shows an abnormal value. Such a point is created, for example, by cosmic rays.



Remove Spike Noise Dialog

The spike noise is searched within a specified range of pixels (Width of Spike Noise). The data point that exceeds the specified noise level (Threshold of Spike Noise) is recognized as the spike noise, and replace by the local average. The number of the removed spike noise will be shown in the dialog.

Display Tags

Displays Image Tags of the front image.

Generated on Thu Mar 29 16:48:24 2018 for Image Processing Utilities by

 1.6.1

Functions

Fourier Transforms

[IPUFFT\(\)](#)
[IPUIFFT\(\)](#)
[IPURealFFT\(\)](#)
[IPURealIFFT\(\)](#)
[IPUPureFFT\(\)](#)
[IPUPureIFFT\(\)](#)
[IPUFFTX\(\)](#)
[IPUIFFTX\(\)](#)
[IPUFFTY\(\)](#)
[IPUIFFTY\(\)](#)
[IPUFFTZ\(\)](#)
[IPUIFFTZ\(\)](#)
[IPUFFTXY\(\)](#)
[IPUIFFTXY\(\)](#)

Cross Correlations

[IPUAutoCorrelation\(\)](#)
[IPUCrossCorrelation\(\)](#)
[IPUFTdCrossCorrelation\(\)](#)

Image Shifting

[IPUImageShift\(\)](#)

Matrix

[IPUMatrixMultiply\(\)](#)
[IPUMatrixInverse\(\)](#)
[IPULUDecomposition\(\)](#)

Informations

[IPUGetVersion\(\)](#)
[IPUGetVersionString\(\)](#)
[IPUGetMKLVersion\(\)](#)
[IPUGetMKLVersionString\(\)](#)

Debugging

[IPUTraceEnabled\(\)](#)
[IPUTraceInit\(\)](#)
[IPUTraceShowTotal\(\) IPUTrace\(\)](#)

VML

[IPUConj\(\)](#)
[IPUSqr\(\)](#)
[IPUAbs\(\)](#)
[IPUInv\(\)](#)
[IPUSqrt\(\)](#)

`IPUInvSqrt()`
`IPUCbrt()`
`IPUInvCbrt()`
`IPUPow2o3()`
`IPUPow3o2()`
`IPUExp()`
`IPUExpm1()`
`IPULn()`
`IPULog10() IPULog1p()`
`IPUCos()`
`IPUSin()`
`IPUTan()`
`IPUCosh()`
`IPUSinh()`
`IPUTanh() IPUAcos()`
`IPUAsin() IPUAtan()`
`IPUAcosh() IPUAsinh()`
`IPUAtanh()`
`IPUErf()`
`IPUErfc()`
`IPUErfinv()`
`IPUFloor()`
`IPUCeil()`
`IPUTrunc()`
`IPURound()`
`IPUNearbyInt()`
`IPURint() IPUAdd()`
`IPUSub() IPUMul()`
`IPUDiv()`
`IPUPow()`
`IPUHypot() IPUAtan2()`
`IPUMulByConj()`
`IPUAbs()`
`IPUCIS()`
`IPUPowx()`
`IPUSinCos()`
`IPUModf()`

File List

Here is a list of all documented files with brief descriptions:

Cross_Correlations.s	
Debugging.s	
Fourier_Transforms.s	
Image_Shifting.s	
Informations.s	
Matrix.s	
VML.s	

Generated on Thu Mar 29 16:48:24 2018 for Image Processing Utilities by



1.6.1

Cross_Correlations.s File Reference

Functions

Image IPUAutoCorrelation (Image imgSrc)	
Auto Correlation.	
Image IPUCrossCorrelation (Image imgTarg, Image imgRef)	Cross
Correlation.	
Image IPUCrossCorrelation (Image imgTarg, Image imgRef, Number eMode, Number fSize)	
Cross Correlation.	
Image IPUFTdCrossCorrelation (ComplexImage cimTarg, ComplexImage cimRef, Number eMode, Number fSize)	
Cross Correlation from Fourier Transformed Images.	

Detailed Description

Function Documentation

Image **IPUAutoCorrelation** (Image *imgSrc*)

Auto Correlation.

Parameters:

imgSrc Source image

Returns:

Auto correlation image

Same as AutoCorrelation()

Image **IPUCrossCorrelation** (Image *imgTarg*, Image *imgRef*)

Cross Correlation.

Parameters:

imgTarg Target image

imgRef Reference image

Returns:

Cross correlation image

Same as CrossCorrelation()

Image **IPUCrossCorrelation** (Image *imgTarg*, Image *imgRef*, Number *eMode*,

Number fSize

)

Cross Correlation.

Parameters:

imgTarg Target image
imgRef Reference image
eMode Cross correlation mode
 1..Conventional, 2..Square Root Moduli, 3..Phase
fSize Mask size to remove high frequency noise.
 fSize is a radius given in pixels in Fourier space, or 0 to 0.5 in terms
 of Nyquist frequency.

Returns:

Cross correlation image

**Image IPUFTdCrossCorrelation (ComplexImage cimTarg,
 ComplexImage cimRef,
 Number eMode,
 Number fSize
)**

Cross Correlation from Fourier Transformed Images.

Parameters:

cimTarg Fourier transformed target image
cimRef Fourier transformed reference image
eMode Cross correlation mode
 1..Conventional, 2..Square Root Moduli, 3..Phase
fSize Mask size to remove high frequency noise.
 fSize is a radius given in pixels in Fourier space, or 0 to 0.5 in terms
 of Nyquist frequency.

Returns:

Cross correlation image

Debugging.s File Reference

Functions

void	IPUTraceInit ()	Initialize Trace Routine.
void	IPUTraceShowTotal ()	Show Total Execution.
void	IPUTraceEnabled (Number isEnabled)	Set Trace Status.
Object	IPUTrace (String strName)	Enter a Trace Block.
Object	IPUTrace (String strName, String strArgs)	Enter a Trace Block with Optional Message.

Detailed Description

These routines show call stacks and an execution profile.

Example:

```
Number MyFunc (Number n)
{
    Object tracer = IPUTrace("MyFunc", ""+n);
    if ( n > 1 ) {
        return n + MyFunc(n - 1);
    } else
    {
        return
n;
    }
}

{
    IPUTraceEnabled(1);
    IPUTraceInit();

    Result(MyFunc(3)+"\n");

    IPUTraceShowTotal();
}
```

This will show the following results:

```
MyFunc(3) {
    MyFunc(2) {
        MyFunc(1) {
            } 0.346692 ms
        } 4.86766 ms
    } 7.39452 ms
6

Total Execution:
    MyFunc : 3 calls, 12.6 ms
```

Function Documentation

void IPUTraceInit ()

Initialize Trace Routine.

example

You should call this before starting a trace.

void IPUTraceShowTotal ()

Show Total Execution.

example

This shows total execution information, if the trace is enabled.

void IPUTraceEnabled (Number isEnabled)

Set Trace Status.

Parameters:

isEnabled 1/0 for enable/disable.

example

Object IPUTrace (String strName)

Enter a Trace Block.

Parameters:

strName Block name to be shown

Returns:

Tracer

example

Object IPUTrace (String strName, String strArgs)

Enter a Trace Block with Optional Message.

Parameters:

strName Block name to be shown

strArgs Any message to be shown

Returns:

Tracer

example

Generated on Thu Mar 29 16:48:24 2018 for Image Processing Utilities by  1.6.1

Fourier_Transforms.s File Reference

Functions

ComplexImage	IPUFFT (ComplexImage cimSrc) Fourier Transform (Origin at Center).
ComplexImage	IPUFFT (ComplexImage cimSrc, ComplexImage cimDst) Fourier Transform (Origin at Center).
ComplexImage	IPUFFT (ComplexImage cimSrc, Number fScale, ComplexImage cimDst) Fourier Transform (Origin at Center).
ComplexImage	IPUIFFT (ComplexImage cimSrc) Inverse Fourier Transform (Origin at Center).
ComplexImage	IPUIFFT (ComplexImage cimSrc, ComplexImage cimDst) Inverse Fourier Transform (Origin at Center).
ComplexImage	IPUIFFT (ComplexImage cimSrc, Number fScale, ComplexImage cimDst) Inverse Fourier Transform (Origin at Center).
ComplexImage	IPURealFFT (Image imgSrc) Fourier Transform with ReallImage.
ComplexImage	IPURealFFT (Image imgSrc, ComplexImage cimDst) Fourier Transform with ReallImage.
ComplexImage	IPURealFFT (Image imgSrc, Number fScale, ComplexImage cimDst) Fourier Transform with ReallImage.
Image	IPURealIFFT (ComplexImage cimSrc) Inverse Fourier Transform to ReallImage.
Image	IPURealIFFT (ComplexImage cimSrc, Image imgDst) Inverse Fourier Transform to ReallImage.
Image	IPURealIFFT (ComplexImage cimSrc, Number fScale, Image imgDst) Inverse Fourier Transform to ReallImage.
ComplexImage	IPUPureFFT (ComplexImage cimSrc, Number fScale, ComplexImage cimDst) Fourier Transform (Origin at Top-Left).
ComplexImage	IPUPureIFFT (ComplexImage cimSrc, Number fScale, ComplexImage cimDst) Inverse Fourier Transform (Origin at Top-Left).
ComplexImage	IPUFFTX (ComplexImage cimSrc, Number fScale, ComplexImage cimDst) 1D Fourier Transform at X axis (Origin at Center)
ComplexImage	IPUFFTX (ComplexImage cimSrc, ComplexImage cimDst) 1D Fourier Transform at X axis (Origin at Center)
ComplexImage	IPUFFTX (ComplexImage cimSrc) 1D Fourier Transform at X axis (Origin at Center)
ComplexImage	IPUIFFTX (ComplexImage cimSrc, Number fScale, ComplexImage cimDst) 1D Inverse Fourier Transform at X axis (Origin at Center)
ComplexImage	IPUIFFTX (ComplexImage cimSrc, ComplexImage cimDst) 1D Inverse Fourier Transform at X axis (Origin at Center)
ComplexImage	IPUIFFTX (ComplexImage cimSrc)

	1D Inverse Fourier Transform at X axis (Origin at Center)
ComplexImage	I_PU_FF_TY (ComplexImage cimSrc, Number fScale, ComplexImage cimDst)
	1D Inverse Fourier Transform at Y axis (Origin at Center)
ComplexImage	I_PU_FF_TY (ComplexImage cimSrc, ComplexImage cimDst)
	1D Fourier Transform at Y axis (Origin at Center)
ComplexImage	I_PU_FF_TY (ComplexImage cimSrc)
	1D Fourier Transform at Y axis (Origin at Center)
ComplexImage	I_PU_IFF_TY (ComplexImage cimSrc, Number fScale, ComplexImage cimDst)
	1D Inverse Fourier Transform at Y axis (Origin at Center)
ComplexImage	I_PU_IFF_TY (ComplexImage cimSrc, ComplexImage cimDst)
	1D Inverse Fourier Transform at Y axis (Origin at Center)
ComplexImage	I_PU_IFF_TY (ComplexImage cimSrc)
	1D Inverse Fourier Transform at Y axis (Origin at Center)
ComplexImage	I_PU_FF_TZ (ComplexImage cimSrc, Number fScale, ComplexImage cimDst)
	1D Fourier Transform at Z axis (Origin at Center)
ComplexImage	I_PU_FF_TZ (ComplexImage cimSrc, ComplexImage cimDst)
	1D Fourier Transform at Z axis (Origin at Center)
ComplexImage	I_PU_FF_TZ (ComplexImage cimSrc)
	1D Fourier Transform at Z axis (Origin at Center)
ComplexImage	I_PU_IFF_TZ (ComplexImage cimSrc, Number fScale, ComplexImage cimDst)
	1D Inverse Fourier Transform at Z axis (Origin at Center)
ComplexImage	I_PU_IFF_TZ (ComplexImage cimSrc, ComplexImage cimDst)
	1D Inverse Fourier Transform at Z axis (Origin at Center)
ComplexImage	I_PU_IFF_TZ (ComplexImage cimSrc)
	1D Inverse Fourier Transform at Z axis (Origin at Center)
ComplexImage	I_PU_FF_TXY (ComplexImage cimSrc, Number fScale, ComplexImage cimDst)
	2D Fourier Transform at XY Plane (Origin at Center)
ComplexImage	I_PU_FF_TXY (ComplexImage cimSrc, ComplexImage cimDst)
	2D Fourier Transform at XY Plane (Origin at Center)
ComplexImage	I_PU_FF_TXY (ComplexImage cimSrc)
	2D Fourier Transform at XY Plane (Origin at Center)
ComplexImage	I_PU_IFF_TXY (ComplexImage cimSrc, Number fScale, ComplexImage cimDst)
	2D Inverse Fourier Transform at XY Plane (Origin at Center)
ComplexImage	I_PU_IFF_TXY (ComplexImage cimSrc, ComplexImage cimDst)
	2D Inverse Fourier Transform at XY Plane (Origin at Center)
ComplexImage	I_PU_IFF_TXY (ComplexImage cimSrc)
	2D Inverse Fourier Transform at XY Plane (Origin at Center)

Detailed Description

Function Documentation

ComplexImage IPUFFT (ComplexImage *cimSrc*)

Fourier Transform (Origin at Center).

Parameters:

cimSrc Source image

Returns:

Fourier transformed image

Same as FFT().

ComplexImage IPUFFT (ComplexImage *cimSrc*,

ComplexImage *cimDst*

)

Fourier Transform (Origin at Center).

Parameters:

cimSrc Source image

cimDst Destination image (can be NULL, or *cimSrc* for in-place calculation)

Returns:

cimDst

If *cimDst* is allocated, this function put a result to it. Otherwise, this function allocates a new destination area and returns the result on it.

ComplexImage IPUFFT (ComplexImage *cimSrc*,

Number *fScale*,

ComplexImage *cimDst*

)

Fourier Transform (Origin at Center).

Parameters:

cimSrc Source image

fScale Scaling factor

cimDst Destination image (can be NULL, or *cimSrc* for in-place calculation)

Returns:

cimDst

Behavior:

```
cimDst = FFT(cimSrc) * fScale
```

Normally, `fScale` is 1.0.

If `cimDst` is allocated, this function put a result to it. Otherwise, this function allocates a new destination area and returns the result on it.

ComplexImage IPUIFFT (ComplexImage `cimSrc`)

Inverse Fourier Transform (Origin at Center).

Parameters:

`cimSrc` Source image

Returns:

Inverse fourier transformed image

Same as `IFFT()`.

ComplexImage IPUIFFT (ComplexImage `cimSrc`, ComplexImage `cimDst`)

Inverse Fourier Transform (Origin at Center).

Parameters:

`cimSrc` Source image

`cimDst` Destination image (can be NULL, or `cimSrc` for in-place calculation)

Returns:

`cimDst`

If `cimDst` is allocated, this function put a result to it. Otherwise, this function allocates a new destination area and returns the result on it.

ComplexImage IPUIFFT (ComplexImage `cimSrc`, Number `fScale`, ComplexImage `cimDst`)

Inverse Fourier Transform (Origin at Center).

Parameters:

`cimSrc` Source image

`fScale` Scaling factor

`cimDst` Destination image (can be NULL, or `cimSrc` for in-place calculation)

Returns:

`cimDst`

Behavior:

```
cimDst = IFFT(cimSrc) * fScale
```

Normally, *fScale* is 1.0/(Number of pixels).

If *cimDst* is allocated, this function put a result to it. Otherwise, this function allocates a new destination area and returns the result on it.

ComplexImage IPURealFFT (Image *imgSrc*)

Fourier Transform with ReallImage.

Parameters:

imgSrc Source image

Returns:

Fourier transformed image

Same as RealFFT().

ComplexImage IPURealFFT (Image *imgSrc*, ComplexImage *cimDst*)

Fourier Transform with ReallImage.

Parameters:

imgSrc Source image

cimDst Destination image (can be NULL)

Returns:

cimDst

If *cimDst* is allocated, this function put a result to it. Otherwise, this function allocates a new destination area and returns the result on it.

ComplexImage IPURealFFT (Image *imgSrc*, Number *fScale*, ComplexImage *cimDst*)

Fourier Transform with ReallImage.

Parameters:

imgSrc Source image

fScale Scaling factor

cimDst Destination image (can be NULL)

Returns:

cimDst

Warning:

Image type of *imgSrc* should be Real 4 or Real 8.

Behavior:

```
cimDst = FFT(imgSrc) * fScale
```

Normally, *fScale* is 1.0.

If *cimDst* is allocated, this function put a result to it. Otherwise, this function allocates a new destination area and returns the result on it.

Image IPURealIFFT (ComplexImage *cimSrc*)

Inverse Fourier Transform to ReallImage.

Parameters:

cimSrc Source image

Returns:

Inverse fourier transformed image

Same as RealIFFT().

Image IPURealIFFT (ComplexImage *cimSrc*, Image *imgDst*)

Inverse Fourier Transform to ReallImage.

Parameters:

cimSrc Source image

imgDst Destination image (can be NULL)

Returns:

cimDst

Warning:

Image type of *imgDst* should be Real 4 or Real 8.

If *imgDst* is allocated, this function put a result to it. Otherwise, this function allocates a new destination area and returns the result on it.

Image IPURealIFFT (ComplexImage *cimSrc*, Number *fScale*, Image *imgDst*)

Inverse Fourier Transform to ReallImage.

Parameters:

cimSrc Source image

fScale Scaling factor

imgDst Destination image (can be NULL)

Returns:

`cimDst`

Warning:

Image type of `imgDst` should be Real 4 or Real 8.

Behavior:

```
cimDst = IFFT(imgSrc) * fScale
```

Normally, `fScale` is $1.0 / (\text{Number of pixels})$.

If `imgDst` is allocated, this function put a result to it. Otherwise, this function allocates a new destination area and returns the result on it.

```
ComplexImage IPUPureFFT ( ComplexImage cimSrc,
                           Number      fScale,
                           ComplexImage cimDst
                         )
```

Fourier Transform (Origin at Top-Left).

Parameters:

`cimSrc` Source image

`fScale` Scaling factor

`cimDst` Destination image (can be NULL, or `cimSrc` for in-place calculation)

Returns:

`cimDst`

Behavior:

```
cimDst = FFT(cimSrc) * fScale
```

Normally, `fScale` is 1.0.

If `cimDst` is allocated, this function put a result to it. Otherwise, this function allocates a new destination area and returns the result on it.

```
ComplexImage IPUPureIFFT ( ComplexImage cimSrc,
                            Number      fScale,
                            ComplexImage cimDst
                          )
```

Inverse Fourier Transform (Origin at Top-Left).

Parameters:

`cimSrc` Source image

`fScale` Scaling factor

`cimDst` Destination image (can be NULL, or `cimSrc` for in-place calculation)

Returns:

`cimDst`

Behavior:

```
cimDst = IFFT(cimSrc) * fScale
```

Normally, `fScale` is $1.0 / (\text{Number of pixels})$.

If `cimDst` is allocated, this function put a result to it. Otherwise, this function allocates a new destination area and returns the result on it.

```
ComplexImage IPUFFTX ( ComplexImage cimSrc,  
                        Number          fScale,  
                        ComplexImage  cimDst  
                      )
```

1D Fourier Transform at X axis (Origin at Center)

Parameters:

`cimSrc` Source image

`fScale` Scaling factor

`cimDst` Destination image (can be NULL, or `cimSrc` for in-place calculation)

Returns:

`cimDst`

Normally, `fScale` is 1.0.

If `cimDst` is allocated, this function put a result to it. Otherwise, this function allocates a new destination area and returns the result on it.

```
ComplexImage IPUFFTX ( ComplexImage cimSrc,  
                        ComplexImage  cimDst  
                      )
```

1D Fourier Transform at X axis (Origin at Center)

Parameters:

`cimSrc` Source image

`cimDst` Destination image (can be NULL, or `cimSrc` for in-place calculation)

Returns:

`cimDst`

If `cimDst` is allocated, this function put a result to it. Otherwise, this function allocates a new destination area and returns the result on it.

```
ComplexImage IPUFFTX ( ComplexImage cimSrc )
```

1D Fourier Transform at X axis (Origin at Center)

Parameters:

`cimSrc` Source image

Returns:

Fourier transformed image

```
ComplexImage IPUIFFTX ( ComplexImage cimSrc,  
                         Number      fScale,  
                         ComplexImage cimDst  
                         )
```

1D Inverse Fourier Transform at X axis (Origin at Center)

Parameters:

cimSrc Source image

fScale Scaling factor

cimDst Destination image (can be NULL, or *cimSrc* for in-place calculation)

Returns:

cimDst

Normally, *fScale* is 1.0/(Number of pixels).

If *cimDst* is allocated, this function put a result to it. Otherwise, this function allocates a new destination area and returns the result on it.

```
ComplexImage IPUIFFTX ( ComplexImage cimSrc,  
                         ComplexImage cimDst  
                         )
```

1D Inverse Fourier Transform at X axis (Origin at Center)

Parameters:

cimSrc Source image

cimDst Destination image (can be NULL, or *cimSrc* for in-place calculation)

Returns:

cimDst

If *cimDst* is allocated, this function put a result to it. Otherwise, this function allocates a new destination area and returns the result on it.

```
ComplexImage IPUIFFTX ( ComplexImage cimSrc )
```

1D Inverse Fourier Transform at X axis (Origin at Center)

Parameters:

cimSrc Source image

Returns:

Inverse Fourier transformed image

```
ComplexImage IPUFFTY ( ComplexImage cimSrc,
```

```
Number      fScale,  
ComplexImage cimDst  
)
```

1D Inverse Fourier Transform at Y axis (Origin at Center)

Parameters:

cimSrc Source image
fScale Scaling factor
cimDst Destination image (can be NULL, or *cimSrc* for in-place calculation)

Returns:

cimDst

Normally, *fScale* is 1.0.

If *cimDst* is allocated, this function put a result to it. Otherwise, this function allocates a new destination area and returns the result on it.

```
ComplexImage IPUFFTY ( ComplexImage cimSrc,  
                        ComplexImage cimDst  
)
```

1D Fourier Transform at Y axis (Origin at Center)

Parameters:

cimSrc Source image
cimDst Destination image (can be NULL, or *cimSrc* for in-place calculation)

Returns:

cimDst

If *cimDst* is allocated, this function put a result to it. Otherwise, this function allocates a new destination area and returns the result on it.

```
ComplexImage IPUFFTY ( ComplexImage cimSrc )
```

1D Fourier Transform at Y axis (Origin at Center)

Parameters:

cimSrc Source image

Returns:

Fourier transformed image

```
ComplexImage IPUIFFTY ( ComplexImage cimSrc,  
                           Number      fScale,  
                           ComplexImage cimDst  
)
```

1D Inverse Fourier Transform at Y axis (Origin at Center)

Parameters:

cimSrc Source image
fScale Scaling factor
cimDst Destination image (can be NULL, or *cimSrc* for in-place calculation)

Returns:

cimDst

Normally, *fScale* is 1.0/(Number of pixels).

If *cimDst* is allocated, this function put a result to it. Otherwise, this function allocates a new destination area and returns the result on it.

**ComplexImage IPUIFFTY (ComplexImage *cimSrc*,
ComplexImage *cimDst*
)**

1D Inverse Fourier Transform at Y axis (Origin at Center)

Parameters:

cimSrc Source image
cimDst Destination image (can be NULL, or *cimSrc* for in-place calculation)

Returns:

cimDst

If *cimDst* is allocated, this function put a result to it. Otherwise, this function allocates a new destination area and returns the result on it.

ComplexImage IPUIFFTY (ComplexImage *cimSrc*)

1D Inverse Fourier Transform at Y axis (Origin at Center)

Parameters:

cimSrc Source image

Returns:

Inverse Fourier transformed image

**ComplexImage IPUFFTZ (ComplexImage *cimSrc*,
Number *fScale*,
ComplexImage *cimDst*
)**

1D Fourier Transform at Z axis (Origin at Center)

Parameters:

cimSrc Source image
fScale Scaling factor
cimDst Destination image (can be NULL, or *cimSrc* for in-place calculation)

Returns:

cimDst

Normally, *fScale* is 1.0.

If *cimDst* is allocated, this function put a result to it. Otherwise, this function allocates a new destination area and returns the result on it.

**ComplexImage IPUFFTZ (ComplexImage *cimSrc*,
ComplexImage *cimDst*
)**

1D Fourier Transform at Z axis (Origin at Center)

Parameters:

cimSrc Source image

cimDst Destination image (can be NULL, or *cimSrc* for in-place calculation)

Returns:

cimDst

If *cimDst* is allocated, this function put a result to it. Otherwise, this function allocates a new destination area and returns the result on it.

ComplexImage IPUFFTZ (ComplexImage *cimSrc*)

1D Fourier Transform at Z axis (Origin at Center)

Parameters:

cimSrc Source image

Returns:

Fourier transformed image

**ComplexImage IPUIFFTZ (ComplexImage *cimSrc*,
Number *fScale*,
ComplexImage *cimDst*
)**

1D Inverse Fourier Transform at Z axis (Origin at Center)

Parameters:

cimSrc Source image

fScale Scaling factor

cimDst Destination image (can be NULL, or *cimSrc* for in-place calculation)

Returns:

cimDst

Normally, *fScale* is 1.0/(Number of pixels).

If *cimDst* is allocated, this function put a result to it. Otherwise, this function allocates a new destination area and returns the result on it.

```
ComplexImage IPUIFFTZ ( ComplexImage cimSrc,
                  ComplexImage cimDst
)
```

1D Inverse Fourier Transform at Z axis (Origin at Center)

Parameters:

cimSrc Source image

cimDst Destination image (can be NULL, or *cimSrc* for in-place calculation)

Returns:

cimDst

If *cimDst* is allocated, this function put a result to it. Otherwise, this function allocates a new destination area and returns the result on it.

```
ComplexImage IPUIFFTZ ( ComplexImage cimSrc )
```

1D Inverse Fourier Transform at Z axis (Origin at Center)

Parameters:

cimSrc Source image

Returns:

Inverse Fourier transformed image

```
ComplexImage IPUFFTXY ( ComplexImage cimSrc,
                  Number            fScale,
                  ComplexImage cimDst
                 )
```

2D Fourier Transform at XY Plane (Origin at Center)

Parameters:

cimSrc Source image

fScale Scaling factor

cimDst Destination image (can be NULL, or *cimSrc* for in-place calculation)

Returns:

cimDst

Normally, *fScale* is 1.0.

If *cimDst* is allocated, this function put a result to it. Otherwise, this function allocates a new destination area and returns the result on it.

```
ComplexImage IPUFFTXY ( ComplexImage cimSrc,
                  ComplexImage cimDst
                 )
```

2D Fourier Transform at XY Plane (Origin at Center)

Parameters:

cimSrc Source image

cimDst Destination image (can be NULL, or *cimSrc* for in-place calculation)

Returns:

cimDst

If *cimDst* is allocated, this function put a result to it. Otherwise, this function allocates a new destination area and returns the result on it.

ComplexImage IPUFFTXY (ComplexImage *cimSrc*)

2D Fourier Transform at XY Plane (Origin at Center)

Parameters:

cimSrc Source image

Returns:

Fourier transformed image

ComplexImage IPUIFFTXY (ComplexImage *cimSrc*,

Number *fScale*,

ComplexImage *cimDst*

)

2D Inverse Fourier Transform at XY Plane (Origin at Center)

Parameters:

cimSrc Source image

fScale Scaling factor

cimDst Destination image (can be NULL, or *cimSrc* for in-place calculation)

Returns:

cimDst

Normally, *fScale* is 1.0/(Number of pixels).

If *cimDst* is allocated, this function put a result to it. Otherwise, this function allocates a new destination area and returns the result on it.

ComplexImage IPUIFFTXY (ComplexImage *cimSrc*,

ComplexImage *cimDst*

)

2D Inverse Fourier Transform at XY Plane (Origin at Center)

Parameters:

cimSrc Source image

`cimDst` Destination image (can be NULL, or `cimSrc` for in-place calculation)

Returns:

`cimDst`

If `cimDst` is allocated, this function put a result to it. Otherwise, this function allocates a new destination area and returns the result on it.

ComplexImage IPUIFFTXY(ComplexImage `cimSrc`)

2D Inverse Fourier Transform at XY Plane (Origin at Center)

Parameters:

`cimSrc` Source image

Returns:

Inverse fourier transformed image

Generated on Thu Mar 29 16:48:24 2018 for Image Processing Utilities by  1.6.1

Image_Shifting.s File Reference

Functions

ComplexImage **IPUImageShift** (Image imgSrc, Number fShiftX, Number fShiftY, Image imgDst)
Translate Image using FFT.

Detailed Description

Function Documentation

```
ComplexImage IPUImageShift( Image    imgSrc,
                            Number   fShiftX,
                            Number   fShiftY,
                            Image    imgDst
)
```

Translate Image using FFT.

Parameters:

imgSrc Source image
fShiftX X shift
fShiftY Y shift
imgDst Destination image (can be NULL, or **imgSrc** for in-place calculation)

Returns:

imgDst

This function shifts an image using the shift theorem of Fourier Transform. Thus, you can specify sub-pixel numbers to **fShiftX** and **fShiftY**.

If **imgDst** is allocated, this function put a result to it. Otherwise, this function allocates a new destination area and returns the result on it.

Example:

```
Image img := GetFrontImage().ImageClone();
img.ConvertToFloat();

// translate (X, Y) = (10.5, 3)
IPUImageShift(img, 10.5, 3, img);

img.ShowImage();
```

Informations.s File Reference

Number	IPUGetVersion
	Get Version Number.
String	IPUGetVersionString
	Get Version String.
Number	IPUGetMKLVersion
	Get MKL Version Number.
String	IPUGetMKLVersionString
	Get MKL Version String.

Detailed Description

Function Documentation

Number IPUGetVersion ()

Get Version Number.

Returns:

Version number

It returns 0XXYYZZWW.

XX: Major Version

YY: Minor Version

ZZ: Revision Number

WW: Build Number

String IPUGetVersionString ()

Get Version String.

Returns:

Version string

Number IPUGetMKLVersion ()

Get MKL Version Number.

It returns 0XXYYZZZZ.

XX: Major Version

YY: Minor Version

ZZZZ: Build Number

String IPUGetMKLVersionString ()

Get MKL Version String.

Returns:

Version string

Generated on Thu Mar 29 16:48:24 2018 for Image Processing Utilities by  1.6.1

Matrix.s File Reference

Image **IPUMatrixMultiply** (Image imgA, Image imgB)
Matrix Multiply.

Image **IPUMatrixInverse** (Image imgA)
Matrix Inverse.

Image **IPULUDecomposition** (Image imgA, Image imgB)
LUDecomposition.

Image **IPUMatrixTranspose** (Image imgA)
Matrix Transpose.

Detailed Description

Function Documentation

**Image IPUMatrixMultiply (Image imgA,
 Image imgB
)**

Matrix Multiply.

Parameters:

imgA 2D image
imgB 2D image

Returns:

Matrix Multiplied 2D image

Same as MatrixMultiply()

Image IPUMatrixInverse (Image imgA)

Matrix Inverse.

Parameters:

imgA 2D image

Returns:

Inversed Matrix (2D image)

Same as MatrixInverse()

**Image IPULUDecomposition (Image imgA,
 Image imgB
)**

LU Decomposition.

Parameters:

imgA 2D image
imgB 2D image

Returns:

Solution X (2D image) of $AX = B$

Image IPUMatrixTranspose(Image imgA)

Matrix Transpose.

Parameters:

imgA 2D image

Returns:

Transposed Matrix (2D image)

Same as MatrixTranspose()

VML.s File Reference

Image **IPUConj** (ComplexImage imgSrc, ComplexImage imgDst)

Conjugation of vector elements.

Image **IPUSqr** (Image imgSrc, Image imgDst) Squaring of vector elements.

Image **IPUAbs** (Image imgSrc, Image imgDst) Absolute value of vector elements.

Image **IPUInv** (Image imgSrc, Image imgDst) Inversion of vector elements.

Image **IPUSqrt** (Image imgSrc, Image imgDst)

Square root of vector elements.

Image **IPUInvSqrt** (Image imgSrc, Image imgDst)

Inverse square root of vector elements.

Image **IPUCbrt** (Image imgSrc, Image imgDst)

Cube root of vector elements.

Image **IPUInvCbrt** (Image imgSrc, Image imgDst)

Inverse cube root of vector elements.

Image **IPUPow2o3** (Image imgSrc, Image imgDst)

Each vector element raised to 2/3.

Image **IPUPow3o2** (Image imgSrc, Image imgDst)

Each vector element raised to 3/2.

Image **IPUExp** (Image imgSrc, Image imgDst) Exponential of vector elements.

Image **IPUExpm1** (Image imgSrc, Image imgDst)

Exponential of vector elements decreased by 1.

Image **IPULn** (Image imgSrc, Image imgDst)

Natural logarithm of vector elements.

Image **IPULog10** (Image imgSrc, Image imgDst)

Denary logarithm of vector elements.

Image **IPULog1p** (Image imgSrc, Image imgDst)

Natural logarithm of vector elements that are increased by 1.

Image **IPUCos** (Image imgSrc, Image imgDst)

Cosine of vector elements.

Image **IPUSin** (Image imgSrc, Image imgDst)

Sine of vector elements.

Image **IPUTan** (Image imgSrc, Image imgDst)

Tangent of vector elements.

Image **IPUCosh** (Image imgSrc, Image imgDst)

Hyperbolic cosine of vector elements.

Image IPUSinh (Image imgSrc, Image imgDst)

Hyperbolic sine of vector elements.

Image IPUTanh (Image imgSrc, Image imgDst)

Hyperbolic tangent of vector elements.

Image IPUAcos (Image imgSrc, Image imgDst)

Inverse cosine of vector elements.

Image IPUAsin (Image imgSrc, Image imgDst)

Inverse sine of vector elements.

Image IPUAtan (Image imgSrc, Image imgDst)

Inverse tangent of vector elements.

Image IPUAcosh (Image imgSrc, Image imgDst)

Inverse hyperbolic cosine of vector elements.

Image IPUAsinh (Image imgSrc, Image imgDst)

Inverse hyperbolic sine of vector elements.

Image IPUAtanh (Image imgSrc, Image imgDst)

Inverse hyperbolic tangent of vector elements.

Image IPU Erf (Image imgSrc, Image imgDst)

Error function value of vector elements.

Image IPU Erfc (Image imgSrc, Image imgDst)

Complementary error function value of vector elements.

Image IPU ErfInv (Image imgSrc, Image imgDst)

Inverse error function value of vector elements.

Image IPUFloor (Image imgSrc, Image imgDst)

Rounding towards minus infinity. **Image IPUCeil** (Image imgSrc, Image imgDst)

Rounding towards plus infinity.

Image IPUTrunc (Image imgSrc, Image imgDst)

Rounding towards zero infinity.

Image IPURound (Image imgSrc, Image imgDst)

Rounding to nearest integer.

Image IPUNearbyInt (Image imgSrc, Image imgDst)

Rounding according to current mode.

Image IPURint (Image imgSrc, Image imgDst)

Rounding according to current mode and raising inexact result exception.

Image IPUAdd (Image imgA, Image imgB, Image imgDst)

Addition of vector elements.

Image IPUSub (Image imgA, Image imgB, Image imgDst)

Subtraction of vector elements.

Image IPUMul (Image imgA, Image imgB, Image imgDst)

Multiplication of vector elements.

Image IPUDiv (Image imgA, Image imgB, Image imgDst)

Division of elements of one vector by elements of the second vector.

Image IPUPow (Image imgA, Image imgB, Image imgDst)

Each vector element raised to the specified power.

Image IPUHypot (Image imgA, Image imgB, Image imgDst)

Square root of sum of squares.

Image IPUAtan2 (Image imgA, Image imgB, Image imgDst)

Four-quadrant inverse tangent of elements of two vectors.

Image IPUMulByConj (ComplexImage imgA, ComplexImage imgB, ComplexImage imgDst)

Multiplication of elements of one vector by conjugated elements of the second vector.

Image IPUAbs (ComplexImage imgSrc, Image imgDst)

Absolute value of vector elements.

Image IPU CIS (Image imgSrc, ComplexImage imgDst)

Complex exponent of vector elements(cosine and sine combined to complex value).

Image IPU Powx (Image imgSrc, Number nB, Image imgDst)

vector elements raised to the constant power

void IPUSinCos (Image imgSrc, Image imgSin, Image imgCos)

Sine and cosine of vector elements.

void IPUModf (Image imgSrc, Image imgVal, Image imgRem)

Integer and fraction parts.

Detailed Description

Function Documentation

```
Image IPUConj ( ComplexImage imgSrc,  
                 ComplexImage imgDst  
               )
```

Conjugation of vector elements.

Parameters:

imgSrc Source image
imgDst Destination image (can be NULL)

Returns:

imgDst

```
Image IPUSqr ( Image imgSrc,  
                 Image imgDst  
               )
```

Squaring of vector elements.

Parameters:

imgSrc Source image
imgDst Destination image (can be NULL)

Returns:

imgDst

```
Image IPUAbs ( Image imgSrc,  
                 Image imgDst
```

)

Absolute value of vector elements.

Parameters:

imgSrc Source image
imgDst Destination image (can be NULL)

Returns:

imgDst

**Image IPUIInv (Image imgSrc,
 Image imgDst
)**

Inversion of vector elements.

Parameters:

imgSrc Source image
imgDst Destination image (can be NULL)

Returns:

imgDst

**Image IPUSqrt (Image imgSrc,
 Image imgDst
)**

Square root of vector elements.

Parameters:

imgSrc Source image
imgDst Destination image (can be NULL)

Returns:

imgDst

**Image IPUIInvSqrt (Image imgSrc,
 Image imgDst
)**

Inverse square root of vector elements.

Parameters:

imgSrc Source image
imgDst Destination image (can be NULL)

Returns:

imgDst

```
Image IPUCbrt ( Image imgSrc,  
                  Image imgDst  
                )
```

Cube root of vector elements.

Parameters:

imgSrc Source image
imgDst Destination image (can be NULL)

Returns:

imgDst

```
Image IPUInvCbrt ( Image imgSrc,  
                     Image imgDst  
                   )
```

Inverse cube root of vector elements.

Parameters:

imgSrc Source image
imgDst Destination image (can be NULL)

Returns:

imgDst

```
Image IPUPow2o3 ( Image imgSrc,  
                     Image imgDst  
                   )
```

Each vector element raised to 2/3.

Parameters:

imgSrc Source image
imgDst Destination image (can be NULL)

Returns:

imgDst

```
Image IPUPow3o2 ( Image imgSrc,  
                     Image imgDst  
                   )
```

Each vector element raised to 3/2.

Parameters:

imgSrc Source image
imgDst Destination image (can be NULL)

Returns:
imgDst

**Image IPUExp (Image imgSrc,
 Image imgDst
)**

Exponential of vector elements.

Parameters:
imgSrc Source image
imgDst Destination image (can be NULL)

Returns:
imgDst

**Image IPUExpm1 (Image imgSrc,
 Image imgDst
)**

Exponential of vector elements decreased by 1.

Parameters:
imgSrc Source image
imgDst Destination image (can be NULL)

Returns:
imgDst

**Image IPULn (Image imgSrc,
 Image imgDst
)**

Natural logarithm of vector elements.

Parameters:
imgSrc Source image
imgDst Destination image (can be NULL)

Returns:
imgDst

**Image IPULog10 (Image imgSrc,
 Image imgDst
)**

Denary logarithm of vector elements.

Parameters:

imgSrc Source image
imgDst Destination image (can be NULL)

Returns:

imgDst

Image **IPULog1p** (**Image** *imgSrc*,
 Image *imgDst*
)

Natural logarithm of vector elements that are increased by 1.

Parameters:

imgSrc Source image
imgDst Destination image (can be NULL)

Returns:

imgDst

Image **IPUCos** (**Image** *imgSrc*,
 Image *imgDst*
)

Cosine of vector elements.

Parameters:

imgSrc Source image
imgDst Destination image (can be NULL)

Returns:

imgDst

Image **IPUSin** (**Image** *imgSrc*,
 Image *imgDst*
)

Sine of vector elements.

Parameters:

imgSrc Source image
imgDst Destination image (can be NULL)

Returns:

imgDst

Image **INPUTan** (**Image** *imgSrc*,
 Image *imgDst*

)

Tangent of vector elements.

Parameters:

imgSrc Source image
imgDst Destination image (can be NULL)

Returns:

imgDst

**Image IPUCosh (Image imgSrc,
 Image imgDst
)**

Hyperbolic cosine of vector elements.

Parameters:

imgSrc Source image
imgDst Destination image (can be NULL)

Returns:

imgDst

**Image IPUSinh (Image imgSrc,
 Image imgDst
)**

Hyperbolic sine of vector elements.

Parameters:

imgSrc Source image
imgDst Destination image (can be NULL)

Returns:

imgDst

**Image IPUTanh (Image imgSrc,
 Image imgDst
)**

Hyperbolic tangent of vector elements.

Parameters:

imgSrc Source image
imgDst Destination image (can be NULL)

Returns:

imgDst

```
Image IPUAcos ( Image imgSrc,  
                  Image imgDst  
                )
```

Inverse cosine of vector elements.

Parameters:

imgSrc Source image
imgDst Destination image (can be NULL)

Returns:

imgDst

```
Image IPUAsin ( Image imgSrc,  
                  Image imgDst  
                )
```

Inverse sine of vector elements.

Parameters:

imgSrc Source image
imgDst Destination image (can be NULL)

Returns:

imgDst

```
Image IPUAtan ( Image imgSrc,  
                  Image imgDst  
                )
```

Inverse tangent of vector elements.

Parameters:

imgSrc Source image
imgDst Destination image (can be NULL)

Returns:

imgDst

```
Image IPUAcosh ( Image imgSrc,  
                  Image imgDst  
                )
```

Inverse hyperbolic cosine of vector elements.

Parameters:

imgSrc Source image
imgDst Destination image (can be NULL)

Returns:
imgDst

**Image IPUAsinh (Image imgSrc,
 Image imgDst
)**

Inverse hyperbolic sine of vector elements.

Parameters:
imgSrc Source image
imgDst Destination image (can be NULL)

Returns:
imgDst

**Image IPUAtanh (Image imgSrc,
 Image imgDst
)**

Inverse hyperbolic tangent of vector elements.

Parameters:
imgSrc Source image
imgDst Destination image (can be NULL)

Returns:
imgDst

**Image IPUErf (Image imgSrc,
 Image imgDst
)**

Error function value of vector elements.

Parameters:
imgSrc Source image
imgDst Destination image (can be NULL)

Returns:
imgDst

**Image IPUErfc (Image imgSrc,
 Image imgDst
)**

Complementary error function value of vector elements.

Parameters:

imgSrc Source image
imgDst Destination image (can be NULL)

Returns:

imgDst

```
Image IPUErflInv ( Image imgSrc,  
                     Image imgDst  
                   )
```

Inverse error function value of vector elements.

Parameters:

imgSrc Source image
imgDst Destination image (can be NULL)

Returns:

imgDst

```
Image IPUFloor ( Image imgSrc,  
                  Image imgDst  
                )
```

Rounding towards minus infinity.

Parameters:

imgSrc Source image
imgDst Destination image (can be NULL)

Returns:

imgDst

```
Image IPUCeil ( Image imgSrc,  
                  Image imgDst  
                )
```

Rounding towards plus infinity.

Parameters:

imgSrc Source image
imgDst Destination image (can be NULL)

Returns:

imgDst

```
Image IPUTrunc ( Image imgSrc,  
                  Image imgDst
```

)

Rounding towards zero infinity.

Parameters:

imgSrc Source image
imgDst Destination image (can be NULL)

Returns:

imgDst

Image IPURound (**Image** *imgSrc*,
 Image *imgDst*
)

Rounding to nearest integer.

Parameters:

imgSrc Source image
imgDst Destination image (can be NULL)

Returns:

imgDst

Image IPUNearbyInt (**Image** *imgSrc*,
 Image *imgDst*
)

Rounding according to current mode.

Parameters:

imgSrc Source image
imgDst Destination image (can be NULL)

Returns:

imgDst

Image IPURint (**Image** *imgSrc*,
 Image *imgDst*
)

Rounding according to current mode and raising inexact result exception.

Parameters:

imgSrc Source image
imgDst Destination image (can be NULL)

Returns:

imgDst

```
Image IPUAdd ( Image imgA,  
                  Image imgB,  
                  Image imgDst  
)
```

Addition of vector elements.

Parameters:

imgA an image
imgB an image
imgDst Destination image (can be NULL)

Returns:

imgDst

```
Image IPUSub ( Image imgA,  
                  Image imgB,  
                  Image imgDst  
)
```

Subtraction of vector elements.

Parameters:

imgA an image
imgB an image
imgDst Destination image (can be NULL)

Returns:

imgDst

```
Image IPUMul ( Image imgA,  
                  Image imgB,  
                  Image imgDst  
)
```

Multiplication of vector elements.

Parameters:

imgA an image
imgB an image
imgDst Destination image (can be NULL)

Returns:

imgDst

```
Image IPUDiv ( Image imgA,  
                  Image imgB,
```

```
Image imgDst
```

```
)
```

Division of elements of one vector by elements of the second vector.

Parameters:

imgA an image

imgB an image

imgDst Destination image (can be NULL)

Returns:

imgDst

```
Image IPUPow ( Image imgA,
```

```
    Image imgB,
```

```
    Image imgDst
```

```
)
```

Each vector element raised to the specified power.

Parameters:

imgA an image

imgB an image

imgDst Destination image (can be NULL)

Returns:

imgDst

```
Image IPUHypot ( Image imgA,
```

```
    Image imgB,
```

```
    Image imgDst
```

```
)
```

Square root of sum of squares.

Parameters:

imgA an image

imgB an image

imgDst Destination image (can be NULL)

Returns:

imgDst

```
Image IPUAtan2 ( Image imgA,
```

```
    Image imgB,
```

```
    Image imgDst
```

```
)
```

Four-quadrant inverse tangent of elements of two vectors.

Parameters:

imgA an image
imgB an image
imgDst Destination image (can be NULL)

Returns:

imgDst

```
Image IPUMulByConj ( ComplexImage imgA,  
                      ComplexImage imgB,  
                      ComplexImage imgDst  
                    )
```

Multiplication of elements of one vector by conjugated elements of the second vector.

Parameters:

imgA an image
imgB an image
imgDst Destination image (can be NULL)

Returns:

imgDst

```
Image IPUAbs ( ComplexImage imgSrc,  
                 Image           imgDst  
               )
```

Absolute value of vector elements.

Parameters:

imgSrc Source image
imgDst Destination image (can be NULL)

Returns:

imgDst

```
Image IPUCIS ( Image           imgSrc,  
                 ComplexImage imgDst  
               )
```

Complex exponent of vector elements(cosine and sine combined to complex value).

Parameters:

imgSrc Source image
imgDst Destination image (can be NULL)

Returns:

imgDst

```
Image IPUPowx ( Image imgSrc,
                 Number nB,
                 Image imgDst
               )
```

vector elements raised to the constant power

Parameters:

imgSrc Source image
nb Constant value for power b.
imgDst Destination image (can be NULL)

Returns:

imgDst

```
void IPUSinCos ( Image imgSrc,
                  Image imgSin,
                  Image imgCos
                )
```

Sine and cosine of vector elements.

Parameters:

imgSrc Source image
imgSin Sine values
imgCos Cosine values

```
void IPUModf ( Image imgSrc,
                 Image imgVal,
                 Image imgRem
               )
```

Integer and fraction parts.

Parameters:

imgSrc Source image
imgVal Integer values
imgRem Fraction values

All Functions

i

Here is a list of all documented functions, variables, defines, enums, and typedefs with links to the documentation:

- i -

- IPUAbs() : [VML.s](#)
- IPUAcos() : [VML.s](#)
- IPUAcosh() : [VML.s](#)
- IPUAdd() : [VML.s](#)
- IPUAsin() : [VML.s](#)
- IPUAsinh() : [VML.s](#)
- IPUAtan() : [VML.s](#)
- IPUAtan2() : [VML.s](#) IPUAtanh() : [VML.s](#)
- [Cross_Correlations.s](#)
- IPUCbrt() : [VML.s](#)
- IPUCeil() : [VML.s](#)
- IPUCIS() : [VML.s](#)
- IPUConj() : [VML.s](#)
- IPUCos() : [VML.s](#)
- IPUCosh() : [VML.s](#)
- IPUCrossCorrelation() : [Cross_Correlations.s](#)
- IPUDiv() : [VML.s](#)
- IPUErf() : [VML.s](#)
- IPUErfc() : [VML.s](#)
- IPUErfinv() : [VML.s](#)
- IPUExp() : [VML.s](#)
- IPUExpm1() : [VML.s](#)
- IPUFFT() : [Fourier_Transforms.s](#)
- IPUFFTX() : [Fourier_Transforms.s](#)
- IPUFFTXY() : [Fourier_Transforms.s](#)
- IPUFFT() : [Fourier_Transforms.s](#) IPUFFTZ() : [Fourier_Transforms.s](#)
- IPUFloor() : [VML.s](#)
- IPUFTdCrossCorrelation() : [Cross_Correlations.s](#)
- IPUGetMKLVersion() : [Informations.s](#)
- IPUGetMKLVersionString() : [Informations.s](#)
- IPUGetVersion() : [Informations.s](#)
- IPUGetVersionString() : [Informations.s](#)
- IPUHypot() : [VML.s](#)
- IPUIFFT() : [Fourier_Transforms.s](#)
- IPUIFFTX() : [Fourier_Transforms.s](#)
- IPUIFFTXY() : [Fourier_Transforms.s](#)
- IPUIFFT() : [Fourier_Transforms.s](#) IPUIFFTZ() : [Fourier_Transforms.s](#)
- [Image_Shifting.s](#)
- IPUInv() : [VML.s](#)
- IPUInvCbrt() : [VML.s](#)
- IPUInvSqrt() : [VML.s](#)
- IPULn() : [VML.s](#)

`IPULog10() : VML.s`

- `IPULog1p() : VML.s`
- `IPULUDecomposition() : Matrix.s`
- `IPUMatrixInverse() : Matrix.s`
- `IPUMatrixMultiply() : Matrix.s`
- `IPUMatrixTranspose() : Matrix.s`
- `IPUModf() : VML.s`
- `IPUMul() : VML.s`
- `IPUMulByConj() : VML.s`
- `IPUNearbyInt() : VML.s`
- `IPUPow() : VML.s`
- `IPUPow2o3() : VML.s`
- `IPUPow3o2() : VML.s`
- `IPUPowx() : VML.s`
- `IPUPureFFT() : Fourier_Transforms.s`
- `IPUPureIFFT() : Fourier_Transforms.s`
- `IPURealFFT() : Fourier_Transforms.s`
- `IPURealIFFT() : Fourier_Transforms.s`
- `IPURint() : VML.s`
- `IPURound() : VML.s`
- `IPUSin() : VML.s`
- `IPUSinCos() : VML.s`
- `IPUSinh() : VML.s`
- `IPUSqr() : VML.s`
- `IPUSqrt() : VML.s`
- `IPUSub() : VML.s`
- `IPUTan() : VML.s`
- `IPUTanh() : VML.s`
- `IPUTrace() : Debugging.s`
- `IPUTraceEnabled() : Debugging.s`
- `IPUTraceInit() : Debugging.s`
- `IPUTraceShowTotal() : Debugging.s`
- `IPUTrunc() : VML.s`